ANALIZA BIAŁEK

100014883

# Laboratory of Analysis

**Client:** Forge Peptides
**Order number:** 100014883
**Name:** BPC-157 / TB-500
**Number of samples:** 1 vial

**Batch number:** ------
**Receive date:** 2026-05-25
**Start date:** 2026-06-05
**Enda date:** 2026-06-08

## Sample identification

| Substance | Method of taking a sample | Form |
|---|---|---|
| BPC-157, TB-500 | Sample provided by customer | White powder in glass vial |


## Results

| Substance | Sample name | Content | Purity |
|---|---|---|---|
| BPC-157 | BPC-157 / TB-500 20MG | 10,38mg | 99% |
| TB-500 | BPC-157 / TB-500 20MG | 10,20mg | 99% |

## Research comment

The test sample was applied in the amount of 3µl to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H20 and Acetonitrile with the addition of 0.1% trifuluoroacetic acid, a peaks of the sought substances appeared. Its concentration, when compared to the reference chromatogram shows the table above. Height the measurement error is + - 0,05 mg. The actual state of the sample sent before the test is presented on the photo below

## Photo of the test sample




**Password**

M2E87P9D

**Analiza Białek sp. z o.o.**
**KRS:** 0000966168
**NIP:** 8992922970
**REGON:** 521751875


**Headquarters and Laboratory:**
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

**Lab technician:**
Analysis and research specialist
mgr Arkadiusz Zając

ANALIZA BIAŁEK

100014883

## Chromatogram




## Method

| Time | H2O + 0,1% TFA | ACN + 0,1% TFA |
|---|---|---|
| 0min | 95% | 5% |
| 2min | 95% | 5% |
| 60min | 5% | 95% |
| 62min | 5% | 95% |

**Password**

M2E87P9D

**Analiza Białek sp. z o.o.**
**KRS:** 0000966168
**NIP:** 8992922970
**REGON:** 521751875

**Headquarters and Laboratory:**
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

**Lab technician:**
Analysis and research specialist
mgr Arkadiusz Zając